

GMH3511 SETS



GMH 3511-SET

COMFORTABLE MEASUREMENT

GMH3531 SETS



GMH 3531-SET125

FOR LABORATORY AND FIELD

GMH3551 SETS



GMH 3551-SET100

WITH DATA LOGGER

GMH 3511-SET

Art. no. 605021

pH-/Redox-/Temperature measuring instrument
Device complete with pH electrode GE 114, GF1T-T3-B-BS,
5x GPH4, 5x GPH7, 2x GPF100

GMH 3511-G125

Art. no. 475740

pH-/Redox-/Temperature measuring instrument
Device complete with pH electrode GE 125 (PT1000)

General:

For comfortable measurement of pH value and temperature. Even easier operation with a menu reduced to 5 points in GMH 3511.

Minimum measurement effort with maintenance-free gel electrodes and automatic temperature compensation.

Specifications:

see GHM 3511

Accessories and spare parts:

GE 114-BNC

Art. no. 604701

pH-electrode

GE 125-BNC

Art. no. 600731

waterproof pH-electrode incl Pt1000, BNC plug

GKK 3001

Art. no. 611605

with cut-outs for 1 device of the GMH 3xxx series and accessories for water analysis (395 x 295 x 106 mm)

GMH 3531-SET125

Art. no. 474240

pH-/Redox-/Temperature measuring instrument
Device complete GE 125, 5x GPH4, 5x GPH7, 2x GPF100,
GKK 3001

GMH 3531-SET100

Art. no. 604591

pH-/Redox-/Temperature measuring instrument
Device complete GE 100, GF1T-T3-B-BS, 5x GPH4, 5x GPH7,
2x GPF100, GKK 3001

General:

Functional scope for elevated demands in the laboratory and field.

Minimum measurement effort with maintenance-free gel electrodes and automatic temperature compensation.

Specifications:

see GHM 3531

Accessories and spare parts:

GE 100-BNC

Art. no. 600704

pH-electrode, BNC plug

GE 125-BNC

Art. no. 600731

waterproof pH-electrode incl Pt1000, BNC plug

GKK 3001

Art. no. 611605

with cut-outs for 1 device of the GMH 3xxx series and accessories for water analysis (395 x 295 x 106 mm)

GMH 3551-SET125

Art. no. 474903

pH-/Redox-/Temperature measuring instrument with logger; Device complete GE 125, 5x GPH4, 5x GPH7, 2x GPF100, GKK 3001

GMH 3551-SET100

Art. no. 475742

pH-/Redox-/Temperature measuring instrument with logger; Device complete GE 100, GF1T-T3-B-BS, 5x GPH4, 5x GPH7, 2x GPF100, GKK 3001

General:

Sets with very good equipment with integrated data logger

Specifications:

see GHM 3551

Accessories and spare parts:

GE 100-BNC

Art. no. 600704

pH-electrode, BNC plug

GE 125-BNC

Art. no. 600731

waterproof pH-electrode incl Pt1000, BNC plug

GKK 3001

Art. no. 611605

with cut-outs for 1 device of the GMH 3xxx series and accessories for water analysis (395 x 295 x 106 mm)